Whatman™ 13mm GD/X Syringe Filters

Product Information sheet

Warning

For research use only.

Not recommended or intended for diagnosis of disease in humans or animals.

Do not use internally or externally in humans or animals.

The Whatman 13mm GD/X Syringe Filters are designed to enable easy filtration of hard-to-filter samples greater than 10 mL. They are available in a wide variety of filter choices with a polypropylene housing.

Each unit contains a prefiltration stack of Whatman Multigrade GMF 150 (10:1 µm) and GF/F (0.7 µm). The combination of this pre-filtration stack and filter housing design enables you to filter highly particulate loaded samples.

Disposable filtration devices provide labor saving efficiency while insuring consistent filter performance when compared to hand assembled reusable filter housings.

This document provides general information on the products listed below. The specifications in the Technical Data section are intended to provide a basis for establishing functional use, as well as for setting quality assurance test performance levels.

Whatman 13 mm GD/X Syringe Filters

Whatman 13 mm GD/X Syringe Filters are designed to enable easy filtration of hard-to-filter samples up to 10 mL.

Filter Media	Typical Application
CA	Aqueous and some organic samples.
GMF	Aqueous and/or organic samples; high loading capacity.
NYL	Aqueous and/or organic samples; hydrophilic
PP	Aqueous and organic samples
PES	Aqueous based samples; low protein binding
PTFE	Organic based samples. Hydrophobic membrane
PVDF	Aqueous and/or organic based samples; low protein binding membrane

Cellulose Acetate (CA) Membrane Considerations:

The Cellulose Acetate membrane utilized in these devices are "surfactant free". This minimizes the risk of filtrate contamination from aqueous bases extractables associated with the membrane. Cellulose Acetate is a good membrane for biological solutions where high loading capacity and low protein binding are important concerns.

Glass Microfiber Filter (GMF) Considerations.

The Glass Microfiber Filter media are produced entirely from the finest

grade of borosilicate glass microfiber and contain no binders either as manufacturing aids or as wet strengthening additives. Glass Microfiber Filters are generally resistant to weakening or disruption of the fibrous matrix by inorganic or organic solutions and have broad chemical compatibility.

Nylon Membrane (NYL) Considerations:

Nylon membrane is hydrophilic and is a good choice for aqueous and/or aqueous organic samples. The membrane offers good chemical resistance to most common HPLC solvents, however it has limited resistance to acids, bases, halogenated hydrocarbons, aldehydes and strong oxidizing agents. The most common application is HPLC sample filtration.

Polypropylene (PP) Membrane Considerations:

This specially selected Polypropylene membrane provides high flow with both solvent and aqueous compatibility. This medium provides a pure, single-material-construction filter device with a broad range of solvent/chemical resistance.

Polysulfone (PES) Membrane Considerations:

This carefully selected Polysulfone membrane meets the special filtration requirements of durability, high temperature resistance, good chemical compatibility, and low protein adsorption. It is particularly suitable for filtration of serum, plasma and tissue culture solutions as well as other protein containing solutions where minimal adsorptive protein loss is desired.

Polytetrafluoroethylene (PTFE) Membrane Considerations:


Polytetrafluoroethylene membrane is hydrophobic and will not allow water to pass without high pressures. Aqueous solutions may be filtered if the membrane is initially "wetted" with alcohol or another appropriate solvent. Polytetrafluoroethylene membrane will stop aqueous aerosols in gas streams.

Polyvinylidene Fluoride (PVDF) Membrane Considerations:

Polyvinylidene Fluoride membrane is a suitable choice for most HPLC Sample Prep applications. The membrane is slightly hydrophobic with low water breakthrough values. It offers good chemical resistance to common HPLC solvents, has low protein binding and negligible extractables.

Technical Data

Whatman 13 mm GD/X Syringe Filters

Dimensions:	21.7 mm x 29.7 mm
Weight:	Approximately 3 grams
Filtration Area:	1.3 cm ²
Maximum Pressure:	100 psi
Housing:	Polypropylene
Volume "Hold Up"	Full housing 0.5 mL with air purge 50 µL (approx)
Flow Direction:	Flow should enter from the inlet

Connectors:	Inlet - Female Luer Lock (FLL) Outlet - Male Slip Luer (MSL)
Autoclaving:	Autoclave at 121°C (131°C max) at 15 psi for 20 minutes (Not recommended for Nylon)
Glass Microfiber:	100% Borosilicate
Prefiltration Media:	GMF 150 10 µm : 1 µm GF/F 0.7 µm
Filter Media:	See chart below

Integrity Test Data

Description	Pore Size Microns Microns	Minimum Bubble Point Psi
Cellulose Acetate	0.2	46
Cellulose Acetate	0.45	29
GF/A	1.6	N/A
GF/B	1.0	N/A
GF/D	2.7	N/A
GF/F	0.7	N/A
Nylon	0.2	46
Nylon	0.45	29
Polypropylene	0.2	15
Polypropylene	0.45	11
Polyethersulfone	0.2	36
Polyethersulfone	0.45	22
Polytetrafluoroethylene*	0.2	15
Polytetrafluoroethylene*	0.45	8.5
Polyvinylidene Fluoride	0.2	45
Polyvinylidene Fluoride	0.45	28

*Bubble point determined with Methanol. All others determined with water.

Chemical Compatibility of Membrane

Solvent	CA	GMF	NYL	PP	PES	PTFE	PVDF
Acetic Acid 5% +	L	R	R	R	R	R	R
Acetic Acid, Glacial	NR	R	L	R	R	R	R
Acetone	NR	R	R	R	NR	R	NR
Acetonitrile	NR	R	R	R	NR	R	R
Ammonia, 6N	+	L	R	R	R	R	L
Amyl Acetate	NR	R	R	R	NR	R	L
Amyl Alcohol	R	R	R	R	R	R	R
Benzene *	R	R	L	L	NR	R	R
Benzyl Alcohol *	L	R	L	R	NR	R	R
Boric Acid	R	R	L	R	R	R	R
Butyl Alcohol	R	R	R	R	R	R	R
Butyl Chloride *	+	R	NR	NR	+	R	R
Carbon Tetrachloride *	NR	R	L	L	NR	R	R
Chloroform *	NR	R	NR	L	NR	R	R
Cyclohexanone	NR	R	NR	R	NR	R	R
Chlorobenzene	+	R		+	+	R	R
Citric Acid	+	R	R	+	+	R	R
Cresol	NR	R	NR	R	NR	R	NR
Cyclohexane	R	R	R	R	R	R	R
Diethyl Acetamide	NR	R	R	R	NR	R	NR
Dimethyl Formamide	NR	R	R	R	NR	R	NR
Dioxane	NR	R	R	R	NR	R	L
DMSO	NR	R	R	R	NR	R	L
Ethanol	R	R	R	R	R	R	R
Ethers	L	R	R	R	R	R	L
Ethyl Acetate	NR	R	R	R	NR	R	L
Ethylene Glycol	L	R	R	R	R	R	R
Formaldehyde	L	R	R	R	R	R	R
Freon TF	R	R	R	R	R	R	R

Solvent	CA	GMF	NYL	PP	PES	PTFE	PVDF
Formic Acid	L	R	NR	R	L	R	R
Hydrochloric Acid (Conc)	NR	R	NR	L	R	R	R
Hydrofluoric Acid	NR	NR	NR	L	+	R	R
Hexane	R	R	R	R	R	R	R
Isobutyl Alcohol	R	R	R	R	R	R	R
Isopropyl Acetate	NR	R	R	R	NR	R	R
Methanol	R	R	R	R	R	R	R
Methyl Ethyl Ketone	L	R	R	R	NR	R	L
Methylene Chloride *	NR	R	NR	L	NR	R	R
Nitric Acid (Conc)	NR	R	NR	NR	NR	R	NR
Nitric Acid, 6N	L	R	NR	L	L	R	L
Nitrobenzene *	NR	R	L	R	L	R	R
Pentane	R	R	R	L	R	R	R
Perchloro Ethylene	R	R	R	R	NR	R	R
Pyridine	NR	R	L	R	NR	R	R
Phenol (0.5%)	R	R	R	R	L	R	R
Sodium Hydroxide, 6N	NR	NR	L	R	R	R	NR
Sulfuric Acid (Conc)	NR	R	NR	NR	NR	R	NR
Tetrahydrofuran	NR	R	R	L	NR	R	R
Toluene *	L	R	L	L	NR	R	R
Trichloroethane *	NR	R	L	R	NR	R	R
Trichloroethylene *	R	R	NR	R	NR	R	R
Water	R	R	R	R	R	R	R
Xylene *	R	R	L	L	NR	R	R

Note: (R = Resistant; LR = Limited Resistance; NR = Non Resistant; + = Insufficient Data; * = Short term resistance of housing).

Operating Instructions

Safety: High pressures can be obtained when using syringes. The smaller the syringe the higher the pressure that can be generated. As a general guide, the following pressures can be obtained by hand with the syringes indicated: 20 mL, 80 psi; 10 mL, 140 psi; 5 mL, 180 psi; 3 mL, 200 psi; 1 mL, 250 psi. Each user should determine the pressure they can generate by hand with a specific size syringe and take appropriate safety precautions not to exceed the recommended rating for the device used. If these limitations are exceeded, bursting of the device may occur resulting in loss of sample or personal injury.

Pre-Filter Media: The Whatman 25mm GD/X Syringe Filters contain a prefiltration stack of Whatman Multigrade GMF -150 (10:1 µm) and Grade GF/F prefilters. Multigrade GMF -150 is a combination of two glass microfiber filters in one. Manufactured from 100% Borosilicate glass, it's construction consists of a coarse layer on top, meshed with a fine layer below. The Grade GF/F Filter will retain fine particles down to 0.7 µm. The combination of pre filters allows you to filter even the most difficult samples with reduced hand force. Compared to an unprotected membrane, the volume of sample filtered can be three to seven times greater.

Efficiency: To maximize filtration throughout, use the largest pore size filter that will provide the required cleanliness. To extend filter life use low flow rates or pressures.

To use with a Syringe:

- 1) Fill the syringe with the solution to be filtered.
- 2) Secure the filled syringe to the FLL inlet of the syringe filter with a twisting motion.
- 3) With the outlet pointed upward, gradually apply thumb pressure to the syringe plunger to initiate flow.
- 4) Continue thumb pressure until all the air in the device is displaced with liquid.
- 5) Once liquid starts to exit syringe filter outlet, stop applying pressure, point the device downward and away from user.
- 6) Position syringe filter over suitable collection container or other

apparatus and apply pressure again to filter sample. Change filters when flow becomes too slow or resistance becomes excessive.

Air Locks: Seriously hamper flow rates. To eliminate, point the outlet of the filter device upward during the initiation of liquid flow.

Bubble Point (BP) Test: Flush the filter device with 1.0 ml or more of the test fluid. After the filter is completely wet and with the outlet pointed upward, apply air under controlled pressure to the inlet until air breaks through the filter and bubbles can be seen exiting the outlet. The pressure at which air passes through the wetted filter is the BP. Refer to table for typical BP values.

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Note: The most popular sizes are listed above. Other sizes are available upon request.

Legal.

GE, imagination at work and GE Monogram are trademarks of General Electric Company.

Whatman is a trademark of GE Healthcare companies.

© 1997–2013 General Electric Company – All rights reserved.

Previously published 1997.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

For your local office contact information, visit

www.gelifesciences.com/contact

GE Healthcare UK Limited
Amersham Place
Little Chalfont, Buckinghamshire,
HP7 9NA, UK

<http://www.gelifesciences.com>

GE Healthcare offices:

GE Healthcare Bio-Sciences AB
Björkgatan 30, 751 84 Uppsala,
Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg,
Germany

GE Healthcare Bio-Sciences Corp.
800 Centennial Avenue, P.O. Box 1327,
Piscataway, NJ 08855-1327,
USA

GE Healthcare Japan Corporation
Sanken Bldg. 3-25-1, Hyakunincho,
Shinjuku-ku, Tokyo 169-0073,
Japan

imagination at work