

**Certified Reference Material**

This certificate is designed in accordance with ISO 17034 and ISO Guide 31. This certified reference material (CRM) was designed, produced and verified in accordance with ISO/IEC 17025, ISO 17034 and a registered quality management system ISO 9001.

**Product Name**

Pentachlorophenol <sup>13</sup>C<sub>6</sub> 100 µg/mL in Cyclohexane

**Product Code**  
DRE-XA15970100CY

**Lot Number**  
H1492820CY

**CAS No.**  
85380-74-1


**Format**  
Solution

**Mol. Weight**  
272.29

**Expiry Date**  
28 Aug 2030

**Mol. Formula**  
<sup>13</sup>C<sub>6</sub>HCl<sub>5</sub>O

**Storage Temp**  
20°C ± 4°C



**CERTIFIED**  
Concentration  
100.01 µg/mL

**CERTIFIED**  
Expanded Uncertainty (U)  
5.01 µg/mL

**Uncertainty**


The certified value(s) and uncertainty(ies) are determined in accordance with ISO 17034 with an 95% confidence level (k=2). Uncertainty is based on the Total Combined Uncertainty, including uncertainties of preparation, purity of neat materials, homogeneity and stability testing. Stability values are based on real evidence opposed to simulation.

The producer certifies that this certified reference material meets the specification stated in this certificate until the expiry date, provided it is stored unopened at the recommended temperature herein. Product warranties for this certified reference material are set out in the terms and conditions of purchase.

| CERTIFIED BY |             | CERTIFIED ON |            |
|--------------|-------------|--------------|------------|
| N. Müller    | 28 Aug 2024 |              | RM Release |



CHROMATOGRAM



| Signal: * DAD1A, Sig=205,2 Ref=off (b) |          |            |        |       |
|--|----------|------------|--------|-------|
| Nr.                                    | RT [min] | Area       | Height | Area% |
| 1                                      | 17       | 2691.42589 | 187.65 | 100   |

**Instrument**

HPLC/DAD

**Column**

ReproSil 100 C18 5 µm 250 x 3 mm

**Method Details**

Eluent A: 100% Acetonitrile

Eluent B: Water + H3PO4

| Time [min] | Eluent A [%] | Eluent B [%] |
|------------|--------------|--------------|
| 0          | 10           | 90           |
| 1          | 10           | 90           |
| 20         | 98           | 2            |
| 25         | 98           | 2            |

**Flow**

1 mL/min

**Method of Preparation**

The certified value is based on gravimetric and volumetric preparation of this CRM. This CRM has been confirmed by the appropriate analytical techniques.

**Batch Information**

Solvent: Cyclohexane, Lot No. 007323CR46, 260.00 mL.

Concentration calculation based on Combined Purity = Chemical Purity \* Isotopic Value.

**Gravimetric Data**

| Compound Name          | Lot No.  | Combined Weight Purity (%) | (mg)   | Chemical Purity (%) | Isotopic Purity (%) |
|------------------------|----------|----------------------------|--------|---------------------|---------------------|
| Pentachlorophenol 13C6 | G1344558 | 97.2                       | 26.751 | 98.8                | 98.4                |

**Intended Use**

This CRM is intended for use in a laboratory as a calibration and quality control standard or in method development for analytical techniques.

**Safety**

Proper precautions should be observed while handling. See Safety Data Sheet.

**Traceability**

The balances used for gravimetric measurements are calibrated with weights traceable to the national standards (DKD). The calibration of

the balances is verified daily internally and annually by an external accredited calibration service. Chromatographic methods are traceable to the International System of Units (SI).

**Homogeneity**

Random replicate samples of the final packaged CRM have been analysed to prove homogeneity compliant with ISO 17034.

**Storage**

The CRM should be stored in the original sealed container at the indicated temperature.

**Instructions for use**

The CRM should be used shortly after opening to avoid concentration changes due to evaporation. It is recommended to use 1 mL as the minimum sample size and if less material is used, to increase the certified uncertainty by a factor of two for half sample and four for a quarter of sample. If storage after opening is necessary, the CRM should be tightly closed and kept from light and moisture. If the CRM was in a sealed ampoule, it should be transferred to a vial with minimum head space. Visit the support section of our website [lgcstandards.com](http://lgcstandards.com) for a series of Dr. Ehrenstorfer Tech Tip videos and frequently asked questions.

LGC Labor GmbH  
Bgm.-Schlosser-Straße 6A  
86199 Augsburg, Germany  
T | +49 821 906080  
F | +49 821 9060888  
E | [dr.ehrenstorfer@lgcgroup.com](mailto:dr.ehrenstorfer@lgcgroup.com)

LGC Labor GmbH is accredited by DAkkS as reference material producer on ISO 17034 (D-RM-19883-01-00).

All characterisation measurements at LGC Labor GmbH are accredited methods to ISO/IEC 17025.

